



Beach Wheelchair

The Melrose beach wheelchair has been developed to assist disabled children access the beach and water.

Large tyres provide a low pressure footprint allowing the Melrose beach wheelchair to float over the softest sand or pebble beach.

Fabricated from aluminium with a high impact nylon mesh for the upholstery, the Melrose beach wheelchair is light enough to toss in the back of the car. The wheels can be easily removed and the back folds down so storage is no problem.

Custom made to fit you and suit the way you push.



BEACH WHEELCHAIR

STANDARD SPECIFICATIONS

- Standard size to suit children
- Lightweight non rust aluminium frame
- Stainless steel bolts and bearings
- Fold down back
- Flip back armrests
- Adjustable footplate bar height
- Easy release wheels and castors - split pin
- Black fabric sling and tension adjustable back made from high impact, quick dry, nylon mesh
- Sun canopy
- 5 Year Frame Warranty



Flip back armrest allows you to get off sideways onto the beach or into the water


Melrose
KIWI CONCEPT CHAIRS

Contact us for a full catalogue of available accessories.
Specifications are subject to change.

January 2014