







Our custom built Shower / Toilet Wheelchair was designed after consultation with clients. It is manufactured from Aluminium and Stainless Steel to last longer in the bathroom environment.

Its very easy to disassemble and the wheels come off making it a

doddle to transport.

A custom made bag can also be ordered to protect it and to keep the parts together when travelling about.





KOTUKU

PORTABLE SHOWER/TOILET WHEELCHAIR

STANDARD SPECIFICATIONS

Aluminium Frame 90805
Seamless padded Self-Skinning Urethane Seat Cushion 17" Toilet Seat can be faced forward, left, right or rear 5" Solid Front Castors with Push Down to Lock Solid Tyres on 20" Composite Wheels
Swing Back, Pull Out Armrests

Silver Sparkle powdercoat finish Swing away, lift out, flip up Foot Rests 120kg Maximum User Weight

OPTIONS AND ACCESSORIES

Rear mounted "Push to Lock" Brakes (these replace the Locking Castors unless otherwise specified)

Carry Bag (made with or without Castors)

Removable Pan Carrier and Pan

Heavier User Weight option

Custom Seat

