

STALLION JR

CHILD'S ADJUSTABLE FIXED FRAME WHEELCHAIR



The Stallion Junior is a fixed frame wheelchair designed for small children.

The Stallion Junior is custom built to your child's specifications for size and rake. It also has adjustable growth features that allow

us to increase the length and width of the frame as they grow.

Built with components and fabrication methods used to build the adult chairs, these wheelchairs can take the rigours of an active growing child.

stallion junior
MELROSE DAY WHEELCHAIR

Melrose
KIWI CONCEPT CHAIRS

STALLION JUNIOR

CHILD'S ADJUSTABLE FIXED FRAME WHEELCHAIR



STANDARD SPECIFICATIONS

Aluminium Frame 90150

Grow-able Width and Length Frame

Adjustable Height and Angle Seat Back

Fixed or Folding Seat Back

Adjustable Height and Angle, Flip-up Foot-plate

Adjustable Seat Height

Adjustable Balance Point

Castors: Roller-blade
4", 5", 6" Solid
4", 5", 6" Wide Solid

Adjustable Camber

Wheels : Quick Release Spoked Wheels on Alloy Rims
16", 18", 20" or
Quick Release Composite Wheels
12", 16", 18", 20"

Tyres: Standard Tread or Smooth

Hard Anodised Push Rims

Push to Lock Brakes

Black Fabric Seat Sling and Back Upholstery

Five Year Frame Warranty

OPTIONS AND ACCESSORIES

Titanium Frame 90150T

Grow-able Width and Length Frame

Castors: 4", 5" Solid Alloy
3", 4", 5" Soft Roll Alloy

Wheels: Composite Wheels
BMX Wheels
Spinergy Wheels
X Core Wheels

Push-rims: Rubber coated - Red or Black
Vinyl
Foam

Spoke Guards: Clear or Black

Removable Stroller Push Handles

Brakes: Composite Scissor Brakes
Ds-Locks Brakes
Hill climber Brakes
Brake Extension Levers

Anti-Tips with Plastic Wheels or Roller-blades

Adjustable Back Upholstery

Adjustable Back Recline (5 position)

Skirt Guards: Fabric Skirt Guards
Fabric covered ABS Skirt Guards
ABS Drop In Skirt Guards

Armrest: Desk or Full Length Pad
T Adjustable Height, with/without Skirt Guards

Fabric or ABS Foot-plate Cover

Cushion, Leg Strap

Impact Guards, Castor Covers

Tray - slides onto Armrests

Bag - fits behind the Seat Back



Stallion Junior swing away footplates option.



Stallion Junior showing the (Black) Growth Bars that allow us to increase the Frame Width. The removable Stroller Handle is also shown attached to the Lower Growth Bar.