

# WOLVERINE

LACROSSE SPORTS WHEELCHAIR



The Melrose Wolverine comes standard in lightweight aluminium, but can also be built in titanium for a super light super responsive frame.

Completely custom built to the athletes specifications for their body size, physical ability and position they play, all our sports

wheelchairs are built by highly experienced fabricators using the highest quality materials and components to maximize their chance of becoming a Lacrosse Super Hero.

Become an "X-Man" in a Melrose Wolverine.



# WOLVERINE

LACROSSE SPORTS WHEELCHAIR

## STANDARD SPECIFICATIONS

Lightweight aluminium frame

Complete custom fit - built to individual specifications for size, rake and camber

Camber 10° - 20°

Quick release 24" spoked wheels on alloy rims

Rollerblade castors

Black fabric sling and back

Cushion

Leg strap and fabric skirt guards

## OPTIONS AND ACCESSORIES

Super lightweight titanium frame

Tapered seat width and back

Ergo seat and back

Double or single swivel anti-tips

Adjustable tension back

Aluminium skirt guards

Rubber coated Hand Rims - Red or Black

Spoke guards - clear or black (to help sticks "tripping" the spokes)

Tyres, smooth, high pressure IRC and Kenda Iron Cap

Light weight 25", 26" or 700c wheels

Super light weight Spinergy wheels - 24", 25", 26" or 700c with colour option carbon cable spokes

Wheel bag to accommodate four 24", 25" or 26" wheels

