

# Bowls Script Form

Fixed or Folding Frame Sports Wheelchair

90875



**To Order** Fill in the frame dimensions and features required. For variations describe in detail on a separate page. If you have any further questions contact us at the address at the bottom of the page.

**SEND FORM**  
If online, click here to send via e-mail


Date _____	Delivery details _____				
Name _____	_____				
Address _____	_____				
_____	_____				
Phone _____	_____				
Fax / Cell _____	Required by _____				
Email _____	(Delivery is normally 4-6 weeks from order confirmation)				
Disability _____	Grading _____	Weight _____	Height _____	Shoe Size _____	T-Shirt Size _____

## MEASUREMENTS (Refer to diagram below)

- |                                  |  |                                       |   |
|----------------------------------|--|---------------------------------------|---|
| <b>A</b> Seat width _____        | <b>F</b> Frame length _____<br><small>(rear of backpost to centre of castor pivot)</small> | <b>K</b> Castor Offset - right _____  | <b>Q</b> <input type="checkbox"/> Fixed Frame <b>90875X</b>               |
| <b>B</b> Seat depth _____        | <b>G</b> Seat to floor - front _____   | <b>M</b> Camber - left (0°-5°) _____  | <b>R</b> <input type="checkbox"/> Folding Frame <b>90875</b>              |
| <b>C</b> Back height _____       | <b>H</b> Seat to floor - rear _____  | <b>N</b> Camber - right (0°-5°) _____ | <b>Back Style</b> <small>(fixed frame wheelchair only)</small>            |
| <b>D</b> Back angle° _____       | <b>I</b> Balance Point _____   | <b>O</b> Wheel offset - left _____    | <input type="checkbox"/> Folding Back <input type="checkbox"/> Fixed Back |
| <b>E</b> Seat to footplate _____ | <b>J</b> Castor Offset - left _____  | <b>P</b> Wheel offset - right _____   |   |

## ADDITIONAL NOTES

(you can type here)




# Bowls Script Form

Fixed or Folding Frame Sports Wheelchair

Name \_\_\_\_\_

### FOOT REST HANGER ANGLE


- 70°
- 80° (std)
- 90°



### FOOT REST WIDTH


Width for Pro T and V style foot rests only

\_\_\_\_\_




**IMPORTANT**  
Foot rest Width  
When measuring across feet you must add 2" (50mm)

### FIXED FRONT FOOT REST




**IMPORTANT**  
The styles in this box are only available on the Fixed Frame wheelchair




- Straight Hangers
- Pro T Hangers
- V Hangers
- Double Bow style foot plate (std)
- D style foot plate (std)
- Angle adjustable foot plate **AD11520**
- 2 piece, flip up plastic foot plates **AD11524**
- 1 piece, flip up angle adjustable alloy foot plate, opens to the left **AD11550**
- 1 piece, flip up angle adjustable alloy foot plate, opens to the right (pictured) **AD11550**
- ABS on foot plate **AD11500**
- 2 piece, flip up angle adjustable alloy foot plates **AD11521**

### DROP IN FOOT REST



- Straight Hangers
- Pro T Hangers
- Double Bow foot plate **AD11580**
- D style foot plate **AD11580**
- Angle adjustable foot plate **AD11520**
- ABS on foot plate **AD11500**

### SWING AWAY FOOT REST



- Straight Hangers
- Pro T Hangers
- V Hangers
- 2 piece, flip up plastic foot plates **AD11560**
- 1 piece, flip up angle adjustable alloy foot plate, opens to the left **AD11555**
- 1 piece, flip up angle adjustable alloy foot plate, opens to the right (pictured) **AD11555**
- 2 piece, flip up angle adjustable alloy foot plates **AD11570**

# Bowls Script Form

Fixed or Folding Frame Sports Wheelchair

Name \_\_\_\_\_

## CASTORS

8" (200 x 50) pneumatic caster with tube wrap around (std)

## FORKS

Standard forks (std)

## PUSH HANDLES

Nil

In back canes (std) (folding frame wheelchair only)

Foldable in back canes **AD14160** (folding or fixed frame wheelchair)

Removable on grab bar **AD14100** (fixed frame wheelchair only)

Dual clamp adj height **AD14120** (fixed frame wheelchair only)

## BRAKES

Nil

Scissor (std)

Superlock **AD14560** (push to lock with brake extensions)

Brake extensions - left **AD14550**

Push to lock (std)

Hill climbers **AD14530**

Brake extensions - right **AD14550**

## ANTI-TIPS

Nil

With plastic wheels **AD15000**

Adjustable anti tips **AD15010**

With roller blades **AD15030**

Kick tube **AD15040**

## ARMRESTS

Nil

T adjustable height with desk pad armrest **AD15260**

ABS skirt guard attached to full length pad armrest **AD15281**

T adjustable height with full length pad armrest **AD15261**

## UPHOLSTERY

Sling (std)

Standard back (std)

No cushion

Bag under chair **AD16600**

Tension adjustable back **AD16010**

Cushion 1" **AD16201**

Netting under chair **AD16650**

Cushion 2" **AD16202**

Bag on back of chair **AD16680**

Cushion 3" **AD16203**

Wheel bag **AD16780**

## SKIRT GUARDS

Nil

ABS drop in **AD16110**

Alloy fixed **AD16130**

Carbon Fibre folding **AD16160**

Fabric skirt guards **AD16100**

ABS fabric covered **AD16120**

Alloy folding **AD16140**

## STRAPS / GLOVES / IMPACT GUARDS

Leg strap (std)

Gloves

Impact guards front hangers **AD16400**

Heel loop straps (std)

(size & style)

Impact guards back bars **AD16420**

## WHEELS

20" BMX (std)

Other: \_\_\_\_\_

## TYRES

20" x 1.75" Smooth tyres with tube wrap around (std)

Other: \_\_\_\_\_

## PUSHRIMS

Anodised (std)

Rubber Coated **AD13850**

Foam Coated **AD13840**

Vinyl Coated **AD13830**

## PUSHRIMS - POSITION

Standard

In close

Onto rim

## ACCESSORIES

Ball holder (tray) **AD17020**

Ball Grabber **AD17021**

# Bowls Script Form

## Folding Frame Sports Wheelchair

Name \_\_\_\_\_

**COLOURS**

- 27 Candy Blue
- 1 Royal Blue
- 72 Glass Blue
- 16 Candy Purple
- 48 Burgundy
- 83 Candy Raspberry
- 88 Fluro Pink
- 28 Candy Red
- 11 Flame
- 57 Super Hot Chilli
- 80 Candy Orange
- 71 Pure Orange
- 68 Crystal Yellow
- 73 Crystal Jade
- 89 Fluro Green
- 70 Candy Emerald
- 2 Hawthorn Green
- 17 Crystal Sapphire

**ATTENTION**

**THIS IS A GUIDE TO COLOURS ONLY, how colours appear on your monitor or when printed out may vary.**

Contact us if you require assistance or a colour sample.

This colour selection is for our aluminum frames only.

We do not recommend powder coating titanium.

Colour availability is subject to change without notice, contact us if you are concerned about the availability of your colour.

- 55 Silver Sparkle
- 12 Silver Orbiter
- 46 Chrome
- 42 Metallic Viper
- 30 Satin Black
- 61 High Gloss Black
- 62 Matt Black
- 63 Texture Matt Black
- 18 Clear
- Bead Blast** (titanium only)
- Custom Colour** (additional cost)
- Description

