

# Condor Bariatric Script

## Bariatric Folding Frame Day Wheelchair

90050B



**To Order** Fill in the frame dimensions and features required. For variations describe in detail on a separate page. If you have any further questions contact us at the address at the bottom of the page.

**SEND FORM**  
If online, click here to send via e-mail

Date \_\_\_\_\_

Name \_\_\_\_\_

Address \_\_\_\_\_

\_\_\_\_\_

Phone \_\_\_\_\_

Fax / Cell \_\_\_\_\_

Email \_\_\_\_\_

Delivery details \_\_\_\_\_

\_\_\_\_\_

Therapist details \_\_\_\_\_

\_\_\_\_\_

Required by \_\_\_\_\_

**TRIAL CHAIR REF:** \_\_\_\_\_

Disability \_\_\_\_\_ Weight \_\_\_\_\_ Height \_\_\_\_\_ Shoe Size \_\_\_\_\_ T-Shirt Size \_\_\_\_\_

### MEASUREMENTS (Refer to diagram)


**Max User Weight 250kg**  
**Cross Arms: 4**

Materials  Aluminium

<b>A</b> Seat width	22" <input type="checkbox"/> 24" <input type="checkbox"/>	<b>D</b> Back angle (0° or 5°)	_____	<b>I</b> Camber (0°, 3°, 5°)	_____
<b>A</b> Custom seat width	_____	<b>E</b> Seat to footplate	_____	<b>J</b> Hanger Angle (70°, 80° or 90°)	_____
<b>B</b> Seat depth	18" <input type="checkbox"/> 20" <input type="checkbox"/> 22" <input type="checkbox"/> 24" <input type="checkbox"/>	<b>F</b> Seat to floor - front	_____	<b>Axle plate</b>	
<b>B</b> Custom seat depth	_____	<b>G</b> Seat to floor - rear	_____	<input type="checkbox"/> Standard <input type="checkbox"/> Amputee	
<b>C</b> Back height	_____	<b>H</b> Balance Point	_____		

### ADDITIONAL NOTES

(you can type here)



# Condor Bariatric Script

Bariatric Folding Frame Day Wheelchair

Name \_\_\_\_\_

**HANGER TYPE**

Swing Away



Drop In




Fixed



**FOOT-PLATE WIDTH**


Width for Pro T and V style foot-plates




When measuring across feet you must add 2" (50mm)

**SWING AWAY FOOTREST**


Straight swing away  
2 flip up




Pro T swing away  
2 flip up




V swing away  
2 flip up




Straight swing away  
1 flip up



Pro T swing away  
1 flip up



V swing away  
1 flip up



Straight Hangers

2 piece flip up, angle adjustable foot-plates (std)

Pro T Hangers

1 piece flip up, angle adjustable, alloy foot-plate, opens to the left **AD11450**

V Hangers

1 piece flip up, angle adjustable, alloy foot-plate, opens to the right (pictured) **AD11450**

**DROP IN FOOTREST**


Straight drop in D  
foot-plate




Pro T drop in D  
foot-plate



Straight drop in  
angle adj



Pro T drop in  
angle adj



Straight Hangers

D style foot-plate (std)


ABS on foot-plate **AD11500**

Pro T Hangers


Angle adjustable foot-plate **AD11520**

**FIXED FOOTREST**


Straight fixed  
2 flip up




Pro T fixed  
2 flip up



V fixed  
2 flip up



Straight fixed  
1 flip up



Pro T fixed  
1 flip up



V fixed  
1 flip up



Straight Hangers

2 piece flip up, angle adjustable foot-plates (std)

Pro T Hangers

1 piece flip up, angle adjustable, alloy foot-plate, opens to the left **AD11450**

V Hangers

1 piece flip up, angle adjustable, alloy foot-plate, opens to the right (pictured) **AD11450**

**ELEVATING LEG REST**

Made for the Swing Away front only

Left **AD11610**

Right **AD11610**

Both **AD11600**

**STUMP SUPPORT**

Made for the Fixed and Swing Away fronts only

Left **AD11650**

Fixed

Right **AD11650**

Removable

Both **AD11651**

# Condor Bariatric Script

## Bariatric Folding Frame Day Wheelchair

Name \_\_\_\_\_

**CASTORS**

- 5" solid alloy (std)     
  6" solid alloy (std)     
  8" floater wide **AD11761**  
 5" x 1.5" alloy **AD11740**     
  6" x 1.5" alloy **AD11742**  
 5" wide soft roll **AD11810**     
  6" wide soft roll **AD11820**

**FORKS**

- Standard forks (std)

**PUSH HANDLES**

- Nil     
  In back canes (std)     
  Stroller D removable **AD14200**  
 Foldable in back canes **AD14160**     
  Dual clamp adj height **AD14120**

**BRAKES**

- Nil     
  Push to lock (long handle)     
  Brake extensions **AD14550**     
  Pull to lock **AD14510**  
 Push to lock (std)     
  Reverse handles     
  Hill climbers **AD14530**     
  Scissor **AD14520**

**ANTI-TIPS**

- Nil     
  Stainless steel with plastic wheels (std)     
  Stainless steel adjustable anti tips **AD15010**  
 Kick tube **AD15040**     
  Stainless steel with rollerblades **AD15030**

**ARMRESTS**

- Nil     
  T adjustable height desk pad armrest with ABS skirt guard **AD15280**  
 T adjustable height with desk pad armrest **AD15260**     
  T adjustable height full pad armrest with ABS skirt guard **AD15281**  
 T adjustable height with full length pad armrest **AD15261**

**UPHOLSTERY**

- Sling (std)     
  Tension adjustable back (std)     
  No cushion     
  Bag under chair **AD16600**  
 No back upholstery     
  Cushion 1" **AD16201**     
  Bag on back of chair **AD16680**  
 Cushion 2" **AD16202**     
  Bag on front left hanger **AD16620**  
 Cushion 3" **AD16203**     
  Bag on front right hanger **AD16620**

**SKIRT GUARDS**

- Nil     
  Fabric skirt guards **AD16100**     
  ABS skirt guards drop in **AD16110**

**STRAPS / BELTS, COVERS / GUARDS**

- Leg strap (std)     
  Fabric foot-plate cover **AD11510**     
  Seat belt - plastic buckle **AD16500**  
 Heel loop straps (std)     
  Castor pivot covers **AD16410**     
  Seat belt - adults car buckle **AD16522**  
 Hanger Impact guards **AD16400**

**TRAY**

- Nil  
 Slides over armrests (reinforced) **AD17001**

# Condor Bariatric Script

Bariatric Folding Frame Day Wheelchair

Name \_\_\_\_\_

**WHEELS****24"** Bariatric Spoked Wheels (std)**TYRES****24"** 24" x 1-3/8" Tread (std) 24" x 1-3/8" Solid (std) 24" x 1-3/8" Kenda HP Tread **AD13710** 24" x 1-3/8" Primo V-Trak Knobbly **AD13715****PUSHRIMS** Anodised (std) Rubber coated (black) **AD13850** Capstan

No. of projections:

Angle:

Length: \_\_\_\_\_

 Anodised (1") **AD13855** Rubber coated (1") **AD13857** Vinyl **AD13830** Foam **AD13840** Vinyl (1") **AD13859****PUSHRIMS - POSITION** Standard In close Onto rim**AXLES** Quick release axles (std)**SPOKE GUARDS** Clear polycarb **AD13900** Black ABS **AD13901**


Name \_\_\_\_\_

**ATTENTION**

**THIS IS A GUIDE TO COLOURS ONLY, how colours appear on your monitor or when printed out may vary.**

Talk to us if you want your Titanium frame coloured or polished. *Colour availability is subject to change without notice.*


**STANDARD & ENABLE CLIENT COLOURS**



**GLASS COLOURS**



**CANDY COLOURS**



**CANDY COLOURS (SPECIAL)**



**FLUORESCENT COLOURS**

