

Limited Edition Script

Fixed Frame Day Wheelchair

90470

Melrose
KIWI CONCEPT CHAIRS

To Order Fill in the frame dimensions and features required. For variations describe in detail on a separate page. If you have any further questions contact us at the address at the bottom of the page.

SEND FORM
If online, click here to send via e-mail

Date _____

Delivery details _____

Name _____

Address _____

Therapist details _____

Phone _____

Fax / Cell _____

Required by _____

Email _____

TRIAL CHAIR REF: _____

Disability _____ Weight _____ Height _____ Shoe Size _____ T-Shirt Size _____

MEASUREMENTS (Refer to diagram below)

Frame Materials

Aluminium (std) Titanium **90470T**

A Seat Width _____

F Frame Length _____

K Castor Pivot Centres _____

Ergo Seat Option **AD10220**

B Seat Length _____

G Seat to Floor - Front _____

L Footplate Width _____

N Ergo Seat - Length _____

C Back Height _____

H Seat to Floor - Rear _____

M Camber (0° - 5°) _____

O Ergo Back - Height _____

D Back Angle° _____

I Balance Point _____

Front Hanger Taper Style

P Ergo Back - Angle° _____


E Seat to Footplate _____

J Castor Pivot Position _____

V Front Pro T Front

ADDITIONAL NOTES

(you can type here)



Limited Edition Script Form

Fixed Frame Day Wheelchair

Name _____

BACK REST

- Folding back (std) Fixed back

FOOTPLATE

- Standard Angle adjustable AD11520 ABS on footplate AD11500 Fabric padded cover AD11510

CASTORS

- Rollerblade (std) 4" solid (std) 5" solid (std)
 4" x 1.25" alloy AD11710 5" solid alloy AD11740
 4" wide soft roll AD11806 5" x 1.5" solid AD11720
 5" x 1.5" alloy AD11740
 5" wide soft roll AD11810

FORKS

- Ultralight bearing forks (std) Frog Legs suspension AD11950 Frog Legs Unitine AD11920

PUSH HANDLES

- Nil Removable AD14100

BRAKES

- Nil Superlock Push to Lock with Brake Extensions AD14560 Scissor (std) Scissor Lock Composite AD14570
 Push to lock (std) Brake extension left AD14550 Ergo Scissor AD14522 D's-Locks AD14600
 Pull to lock AD14510 Brake extension right AD14550 Ergo Scissor w/spring AD14524 Hill climbers AD14530

ANTI-TIPS

- Nil Pair of single anti-tips with plastic wheels AD15000
 Pair of single anti-tips with rollerblades AD15030

ARMRESTS

- Nil ABS skirt guard attached to full length pad armrest AD15281

UPHOLSTERY

- Sling (std) Standard back (std) No cushion Bag under chair AD16600
 Tension adjustable back AD16010 Cushion 1" AD16201 Bag on back of chair AD16680
 Cushion 2" AD16202 Bag on right front hanger AD16620
 Cushion 3" AD16203 Bag on front left hanger AD16620
 Wheel bag AD16780

SKIRT GUARDS

- Nil ABS drop in AD16110 Alloy fixed AD16130 Carbon Fibre fixed AD16161
 Fabric skirt guards AD16100 ABS fabric covered AD16120 Alloy folding AD16140 Carbon Fibre folding AD16160

STRAPS, IMPACT GUARDS

- Leg strap (std) Impact guards front hangers AD16400 Impact guards grab bar AD16420 Impact guards made only if suitable after fitting anti tips, push handles etc.
 Impact guards bottom bar AD16422

Name _____

WHEELS

- 24"**
- Spoked Wheels (std)
 - Spinergy Day Wheels (18 spoke) **AD13300** (Choose your spoke colour) **Bk R Bl W Y**
 - Spinergy LX Wheels (12 spoke) **AD13310** (Choose your spoke colour) **Bk R Bl W Y Gr Pk O**
 - Spinergy LX Wheels with FlexRim (12 spoke - black only) **AD13315**
 - X Core Wheels (5 spoke) **AD13108**
 - Spider Wheels **AD13100**
 - Topolino Wheels **AD13350**
 - Dino Wheels **AD13360** (Choose your spoke colour) **Bk R Bl Sv Go Gr Pk O LB Pu**
 - Glance Wheels **AD13400** (Wheel Type)
 - BMX Wheels (std)

- 25"**
- Spoked Wheels (std)
 - Spinergy Day (18 spoke) **AD13300** (Choose your spoke colour) **Bk R Bl W Y**
 - Spinergy LX (12 spoke) **AD13310** (Choose your spoke colour) **Bk R Bl W Y Gr Pk O**
 - Spinergy LX Wheels with FlexRim (12 spoke - black only) **AD13315**
 - X Core Wheels (3 spoke) **AD13200**
 - Topolino Wheels **AD13350**
 - Dino Wheels **AD13360** (Choose your spoke colour) **Bk R Bl Sv Go Gr Pk O LB Pu**
 - Glance Wheels **AD13400** (Wheel Type)
 - BMX Wheels (std)

- 26"**
- Spoked Wheels (std)
 - Spinergy Day (18 spoke) **AD13300** (Choose your spoke colour) **Bk R Bl W Y**
 - Spinergy LX (12 spoke) **AD13310**
 - BMX Wheels (std)

Spare Wheel(s) Type: _____
Tyres: _____ Pushrims: _____ Size: _____

TYRES

- 24"**
- 24" x 1" HP Smooth (std)
 - 24" x 1" HP Tread (std)
 - 24" x 1" Solid (std)
 - 24" x 1-3/8" Tread (std)
 - 24" x 1-3/8" Solid (std)
 - 24" x 1-3/8" Kenda HP Tread **AD13710**
 - 24" x 1-3/8" Knobbly V-Trak **AD13715**
 - 24" x 1" Kenda Ironcap Blue **AD13700**
 - 24" x 1" Kenda Konstrictor **AD13720**
 - 24" x 1" Schwalbe Right Run **AD13740**
 - 24" x 1" Schwalbe Marathon Evo **AD13750**
 - BMX Tread (std)

- 25"**
- 25" x 1" HP Smooth (std)
 - 25" x 1" Solid (std)
 - 25" x 1" Kenda Ironcap Blue **AD13700**
 - 25" x 1" Kenda Konstrictor **AD13720**
 - 25" x 1" Schwalbe Right Run **AD13740**
 - 25" x 1" Schwalbe Marathon Evo **AD13750**
 - BMX Tread (std)

- 26"**
- 26" x 1" HP Smooth (std)
 - 26" x 1-3/8" Tread (std)
 - 26" x 1" Solid (std)
 - 26" x 1" Schwalbe Right Run **AD13740**
 - 26" x 1" Schwalbe Marathon Evo **AD13750**
 - BMX Tread (std)

AXLES

- Quick release axles (std)
- Quad release **AD12500**

PUSHRIMS

- Anodised (std)
- Anodised (1") **AD13855**
- Foam **AD13840**
- Natural Fit Lite **AD13860**
- Vinyl **AD13830**
- Vinyl (1") **AD13859**
- Surge **AD13864**
- Max Greppe **AD13870**
- Rubber coated (black) **AD13850**
- Rubber coated (1") **AD13857**
- Surge Lite **AD13864**
- Spinergy FlexRim - see Wheels

PUSHRIMS - POSITION



















- Standard
- In close
- Onto rim

SPOKE GUARDS

- Nil
- Clear polycarb **AD13900**
- Black ABS **AD13901**

Name _____











COLOURS

- 27 Candy Blue 
- 1 Royal Blue 
- 72 Glass Blue 
- 16 Candy Purple 
- 48 Burgundy 
- 83 Candy Raspberry 
- 88 Fluro Pink 
- 28 Candy Red 
- 11 Flame 
- 57 Super Hot Chilli 
- 80 Candy Orange 
- 71 Pure Orange 
- 68 Crystal Yellow 
- 73 Crystal Jade 
- 89 Fluro Green 
- 70 Candy Emerald 
- 2 Hawthorn Green 
- 17 Crystal Sapphire 

ATTENTION

THIS IS A GUIDE TO COLOURS ONLY, how colours appear on your monitor or when printed out may vary.

Contact us if you require assistance or a colour sample.
This colour selection is for our aluminum frames only.
We do not recommend powder coating titanium.
Colour availability is subject to change without notice, contact us if you are concerned about the availability of your colour.

- 55 Silver Sparkle 
- 12 Silver Orbiter 
- 46 Chrome 
- 42 Metallic Viper 
- 30 Satin Black 
- 61 High Gloss Black 
- 62 Matt Black 
- 63 Texture Matt Black 
- 18 Clear 
- Bead Blast** (titanium only) 
- Custom Colour** (additional cost)
Description