

AIISLE CHAIR

FOLDABLE AISLE ONBOARD LIFT CHAIR - 90006



The Melrose Folding Lift Aisle Chair is an onboard aisle chair that provides mobility about the cabin during long flights.

It collapses down to 355mm height, 390mm width and 720mm length, this allows it to be safely stowed into most aircraft lockers. The seat lift hoists the user from 470mm up to 700mm allowing easier access over non folding aircraft seat armrests.

The frame is hand made from aircraft grade aluminium and is lightweight and rust resistant. Many of our components are made inhouse to ensure we maintain exact tolerances from the highest quality materials. Flame retardant soft trim materials can be used and can be certified to FAR23.853(a) and FAR25.853(a) App F Part I(a)I(ii) if required.

Aisle Chair

Melrose
KIWI CONCEPT CHAIRS

AISLE CHAIR

FOLDABLE AISLE ONBOARD LIFT CHAIR - 90006

Weight Limit: 150kg

Approx Weight: 21.5kg

STANDARD SPECIFICATIONS

Aluminium Frame 90006

Large range of powder coated colours

Dimensions:

Seat Height from floor: 470mm up to 700mm

Seat Width: 345mm

Seat Length: 430mm

Seat Back Height from cushion: 445mm - *handles extend above this*

Chair width: 390mm

Weight: 21.5kg

Folded Dimensions:

Height: 355mm

Width: 390mm

Length: 720mm

Max user weight: 150kg

Steinco - individual locking 74mm twin castor wheels

Foot operated push to lock brakes

Height and width adjustable footplate

Restraining straps with composite clasp

Five-year frame warranty

OPTIONS AND ACCESSORIES

Seat-belt type buckles on restraining straps

Airline or Airport logo printed on seat back

We can print your logo on the back of the aisle chair.

