

# BOWLS

BOWLS SPORTS WHEELCHAIR



Built to comply with New Zealand Bowls regulations for using wheelchairs on the green. But always get permission from the green-keeper.

Bowls holders can be mounted above or below the chair. The low ride height puts you closer to the green and wider castor

positions allow your feet to clear the castor wheels. Wheel/frame offsets can also be built in to help place you over the bowl.

Built as a folding or fixed frame to aid storage or to provide a robust rigid platform to bowl. What ever you can think of to improve your game, we can build for you.

# BOWLS

## BOWLS SPORTS WHEELCHAIR

### STANDARD SPECIFICATIONS

Lightweight Aluminium Frame

Complete Custom Fit – built to individual specifications for Size, Rake and Camber

Fixed or Folding Frames as standard

Adjustable Back Height in Fixed Frame

Adjustable Footplate Height

Push to Lock Brakes

Quick Release, 20" x 1.75" Wheels

*NZ Bowls Regulations* – Custom 20" x 3" Innova Tyre with 200mm x 75mm Castors

Foot Strap

Black Fabric Sling and Back

5 Year Frame Warranty

### OPTIONS AND ACCESSORIES

Super Lightweight Titanium Frame

Pushrims: Rubber Coated – Red or Black  
Vinyl  
Foam

Push Handles: In Back Canes  
Removable (Only on Fixed Frame Wheelchair)  
Fold Down

Brakes: Hill Climber Brakes  
Scissor Brakes  
Brake Extension Levers

Anti-Tips with Plastic Wheels or Rollerblades

Spoke Guards: Clear or Black

Armrest: Tubular Swing Away  
Desk or Full Length Pad  
T Adjustable Height, with/without Skirt Guards

Adjustable Back Upholstery

Skirt Guards: Fabric Skirt Guards  
Fabric covered ABS Skirt Guards  
ABS Drop In Skirt Guards  
Alloy or Carbon Fibre Skirt Guards

Under Chair Bag or Netting

Ball Holder

Fabric or ABS Footplate Cover

Cushion, Impact Guards, Castor Covers

Tapered Seat Width and Back

Leg Strap, Waist Strap

Wheel Bag to accommodate four 24", 25" or 26" Wheels

