

# GAZELLE 2.0

TENNIS SPORTS WHEELCHAIR



*Adjustable position Tennis Wheelchairs*

*Childs Tennis Wheelchairs*

So you need to get to every ball, you need a stable platform to make shots from, and you need your chair to respond as if it was an extension of your body.

The new Gazelle has a light-weight frame designed to increase first push response. It has four contact points on court for a very stable ride, and because its custom built for you, it fits like a glove.

Built from aluminium or titanium your computer aided designed frame will offer maximum rigidity. This directs all kinetic pushing forces directly to the wheels with no energy lost in frame flexibility.

Custom built and with no adjustment matrix's, the Melrose Gazelle is a simple, uncluttered, very light and durable Tennis Wheelchair.



# GAZELLE

TENNIS SPORTS WHEELCHAIR

## STANDARD SPECIFICATIONS

Lightweight Aluminium 3 or 4 Wheel Frame

Complete custom fit - built to individual specifications for Size, Rake and Camber

Camber 10° - 20°

Quick Release 24", 25" or 26" Spoked Wheels on Alloy Rims

Smooth High Pressure Tyres

Rollerblade Castors

Black Fabric Sling and Back

5 Year Frame Warranty

## OPTIONS AND ACCESSORIES

Super Lightweight Titanium 3 or 4 Wheel Frame

Tapered Seat Width and Back

Custom and angle adjustable footplate options

Double or Single Swivel Anti-Tips

Adjustable Back Upholstery

Skirt Guards: Fabric Skirt Guards  
Alloy or Carbon Fibre Skirt Guards

Straps and Bindings for Hips and Legs

Cushion

Smooth, High Pressure, Asymmetric Tyres

Super light weight Spinergy Wheels - 24", 25" or 26" with colour  
option Carbon Cable Spokes

Wheel Bag to accommodate four 24", 25" or 26" Wheels

