

HAWK ULTRA

FIXED FRAME DAY WHEELCHAIR



The Melrose Hawk ULTRA delivers stunning performance once we place the standard Hawk's Titanium frame on a diet - when we toss out the adjustable seating position matrix - the frame rewards us with a more immediate responsive ride.

Custom built to the individuals specification for seat width, length and height - it has the comfort to match the performance.

When you know that 'less is more'.



HAWK ULTRA

FIXED FRAME DAY WHEELCHAIR

Weight Limit: 100kg max

Approx Weight: 3.9kg - Titanium standard 16" frame without wheels

STANDARD SPECIFICATIONS

Titanium Frame 90420UL

Adjustable Back Height

Seat Width: 13" - 17"

Seat Depth: 13" - 17"

Fixed Back

Fixed Camber Bar Position

Camber Bar: 0° - 5°

Adjustable Foot-plate Height

Ultra-lite Castor Forks

Castors: Roller-blade
4", 5", 6" Solid
4", 5", 6" Wide Solid

Quick release 24" Spoked Wheels on Alloy Rims

Tyres: 24 x 1 High Pressure Tread or Smooth
24 x 1 Solid

Brakes: Push to Lock Brakes
Alloy Scissor Brakes

Black Fabric Sling

Black Fabric Back with Leather Trim

Leg Strap

Three Year Frame Warranty

OPTIONS AND ACCESSORIES

Castors: 4", 5", 6" Solid Alloy
4", 5", 6" Soft Roll
6", 8" Pneumatic

Folding Back

Frog Leg Forks - with standard pivots only

Quad Release Axles

Wheels: Spoked 24", 25", 26"
BMX Wheels
Xcore Wheels
Spinergy Wheels
Loop Wheels
Soft Wheels

Push-rims: Rubber Coated - Red or Black
Vinyl
Foam
CarboLife
Natural Fit

Push handles: In Back Canes
Removable
Fold Down

Brakes: Composite Push to Lock Brakes
Composite Scissor Brakes
Eagle Scissor Brakes
Ds-Locks Brakes
Brake Extension Levers

Anti-Tips with Plastic Wheels or Roller-blades

Spoke Guards: Clear or Black

Adjustable Back Upholstery

Skirt Guards: Fabric Skirt Guards
Fabric covered ABS Skirt Guards
ABS Drop In Skirt Guards
Alloy or Carbon Fibre Skirt Guards

Under Chair Bag

Fabric or ABS Foot-plate Cover

Cushion, Impact Guards, Castor Covers



Castors on our super lightweight fork and bullet pivot system.