

HAWK X

FIXED FRAME DAY WHEELCHAIR



The Melrose Hawk X is our New Generation Hawk and comes standard with an Ultra-Lightweight frame that is compact with no under-frame. The lack of under-frame makes lifting and loading much easier once the rear wheels are popped off.

It is custom built to the individuals specification for seat width

and length for a safe comfortable ride.

Our revised Adjustable Camber Bar position and Castor Pivots allow you to alter your centre of gravity and seat height. This makes it suitable for new wheelchair users that are unsure what seating position will work best for them.



HAWK X

FIXED FRAME DAY WHEELCHAIR

Weight Limit: 100kg max

Approx Weight: 5.5kg - Titanium standard 16" frame without wheels

STANDARD SPECIFICATIONS

Titanium Frame 90420X

Seat Width: 13" - 18"

Adjustable Back Height and Angle

Fixed or Folding Back

Adjustable Foot-plate Height

Adjustable Seat Height

Adjustable Camber Bar Position

Tyres: 24 x 1 High Pressure Tread or Smooth
24 x 1 Solid

Wheels: Quick release 24" Spoked Wheels on Alloy Rims

Castors: Roller-blade
4", 5", 6" Solid
4", 5", 6" Wide Solid

Brakes: Push to Lock Brakes
Alloy Scissor Brakes

Black Fabric Sling

Black Fabric Back with Leather Trim


Leg Strap

Five Year Frame Warranty

A Castor with Adjustable Positioning is necessary when adjusting the Squeeze and Rake.



Standard Adjustable Castor Pivot with Curved Fork.



Standard Adjustable Camber Bar.

OPTIONS AND ACCESSORIES

Castors: 4", 5", 6" Solid Alloy
4", 5", 6" Soft Roll
6", 8" Pneumatic

Forks: Frog Leg Forks
Ultra-lite Forks

Quad Release Axles

Wheels: Spoked 24", 25", 26"
BMX Wheels
Xcore Wheels
Spinergy Wheels

Push-rims: Rubber Coated - Red or Black
Vinyl
Foam
CarboLife
Natural Fit

Brakes: Composite Push to Lock Brakes
Composite Scissor Brakes
Eagle Scissor Brakes
Ds-Locks Brakes
Brake Extension Levers

Push handles: In Back Canes
Removable
Fold Down

Anti-Tips with Plastic Wheels or Roller-blades

Spoke Guards: Clear or Black

Armrest: Tubular Swing Away
Desk or Full Length Pad
T Adjustable Height, with/without Skirt Guards

Adjustable Back Upholstery

Skirt Guards: Fabric Skirt Guards
Fabric covered ABS Skirt Guards
ABS Drop In Skirt Guards
Alloy or Carbon Fibre Skirt Guards

Under Chair Bag

Fabric or ABS Foot-plate Cover

Cushion, Impact Guards, Castor Covers

Camber Bar: 0° - 5°