

KEA

CHILD'S ADJUSTABLE FIXED FRAME WHEELCHAIR



The Melrose Kea is a lightweight aluminium wheelchair that is nice and compact with no under-frame.

Although it is custom built to the child's specifications, we allow room for the wheelchair to grow several inches in length, width

and height. This allows the wheelchair to be used for longer as the child grows.

The adjustable frame also lets you alter the centre of gravity and seat height allowing users to alter their positioning as they get older and more confident.



STANDARD SPECIFICATIONS

Aluminium Frame 90700

Grow-able Width and Length Aluminum Frame

Adjustable Height and Angle Seat Back

Fixed or Folding Seat Back

Adjustable Height and Angle Foot-plate

Adjustable Seat Height

Adjustable Camber Bar Position

Adjustable Balance Point

Push to Lock Brakes

Tyres: Standard Tread or Smooth

Castors: Roller-blade
4", 5" Solid Composite
4", 5" Wide Solid Composite

Wheels : Quick Release Spoked Wheels on Alloy Rims
16", 18", 20" or
Quick Release Composite Wheels
12", 16", 18", 20"

Hard Anodised Push Rims

Black Fabric Seat Sling and Back Upholstery

Five Year Frame Warranty

OPTIONS AND ACCESSORIES

Grow-able Width and Length Titanium Frame

Castors: 4", 5" Solid Alloy
3", 4", 5" Soft Roll Alloy

Wheels: Composite Wheels
BMX Wheels
Spinergy Wheels
X Core Wheels

Push-rims: Rubber coated - Red or Black
Vinyl
Foam

Spoke Guards: Clear or Black

Removable Stroller Push Handles

Brakes: Composite Scissor Brakes
Ds-Locks Brakes
Hill climber Brakes
Brake Extension Levers

Anti-Tips with Plastic Wheels or Roller-blades

Adjustable Back Upholstery

Leg Strap

Armrest: Desk Pad
T Adjustable Height, with/without Skirt Guards

Skirt Guards: Fabric Skirt Guards
Fabric covered ABS Skirt Guards
ABS Drop In Skirt Guards

Fabric or ABS Foot-plate Cover

Cushion

Impact Guards, Castor Covers

Tray - slides onto Armrests

Bag - fits behind the Seat Back



*Recline Option:
Multi position (5) seat
back recliner*



Grow-able options showing the Length and Width adjustments.