

KORIMAKO

FOLDING FRAME DAY WHEELCHAIR



Owning a folding wheelchair needn't exclude you from the stylish and lightweight variations that fixed frame chairs share.

The Korimako is custom built to the clients size and shape allowing us to remove the complexity and weight of the adjustable position parts.

Foot-plates can be made at any angle and can be one piece drop in or swing away flip up. The variations are endless.

Well engineered moving parts ensure the frame reliability and durability of this snazzy folding wheelchair.

KORIMAKO

FOLDING FRAME DAY WHEELCHAIR

Weight Limit:	100kg max, Custom weight option also available
Approx Weight:	8.0kg - Titanium standard 16" frame without wheels 8.9kg - Aluminium standard 16" frame without wheels

STANDARD SPECIFICATIONS

Aluminium Frame 90055

Foot-plates: Swing Away
Drop In
Fixed Flip-Up

Back Canes with Push Handles

Castors: Roller-blade
4", 5", 6" Solid
4", 5", 6" Wide Solid

Quick release 24" Spoked Wheels on Alloy rims

Tyres: 24 x 1 High Pressure Tread or Smooth
24 x 1 Solid

High Mounted Push to Lock Brakes

Black Fabric Sling and Back with Leather trim

Leg Strap or Heel Loops

Three Year Frame Warranty

OPTIONS AND ACCESSORIES

Titanium Frame 90055T

Foot-plate: Specialised Foot-plate options

Seat Backs: Splay Back
Angle Adjustable Back
Folding Back

Castors: 4", 5", 6" Solid Alloy
4", 5", 6" Soft Roll
6", 8" Pneumatic

Frog Leg forks

Quad release Axles

Wheels: Spoked 16", 18", 20", 22", 24", 25", 26"
BMX Wheels
Composite Wheels
Xcore Wheels
Spinergy Wheels
Loop Wheels
Soft Wheels

Push-rims: Rubber coated - Red or Black
Vinyl
Foam
CarboLife
Natural Fit

Spoke Guards: Clear or Black

Push Handles: Adjustable height

Brakes: Composite Scissor Brakes
Eagle Scissor Brakes
Ds-Locks Brakes
Hill Climber Brakes
Brake Extension Levers

Anti-Tips with Plastic Wheels or Roller-blades

Armrest: Tubular Swing Away
Desk or Full Length Pad
T Adjustable height, with/without Skirt Guards

Adjustable Back Upholstery


Skirt Guards: Fabric Skirt Guards
Fabric covered ABS Skirt Guards
ABS Drop In Skirt Guards

Bag under Chair or Back Pack

Cushion

Fabric or ABS Foot-plate Cover

Impact Guards, Castor Covers



Swing Away, Fixed Flip-Up and Drop In Hanger Foot-plate options