

BASKETBALL SPORTS WHEELCHAIR

ORION



Do you need a light-weight super responsive frame designed to increase first push maneuverability?

Built from aluminium or titanium your computer aided designed frame will offer maximum rigidity. This directs all kinetic pushing forces directly to the wheels with no energy lost in frame flexibility.

We custom build each wheelchair to ensure a perfect fit, and to make sure that the athlete gets the very best out of their equipment. We have some of the very best fabricators in the country producing our wheelchairs, and only the best components are used, if they don't meet our high standards then we make them ourselves.



ORION

BASKETBALL SPORTS WHEELCHAIR

STANDARD SPECIFICATIONS

Lightweight Aluminium Frame

Complete custom fit - built to individual specifications for Size, Rake and Camber

Camber 5° - 20°

Quick Release 24", 25" or 26" Spoked Wheels on Alloy Rims

Smooth High Pressure Tyres

Rollerblade Castors

Black Fabric Sling and Back

5 Year Frame Warranty

OPTIONS AND ACCESSORIES

Super Lightweight Titanium Frame

Tapered Seat Width and Back

Double or Single Swivel Anti-Tips

Skirt Guards: Fabric Skirt Guards
Alloy or Carbon Fibre Skirt Guards

Adjustable Back Upholstery

Straps and Bindings for Hips and Legs

Cushion

Smooth, High Pressure, Asymmetric Kenda Kaliente or Schwalbe Speedair Tyres

Super light weight Spinergy Wheels - 24", 25", 26" or 700C with colour option Carbon Cable Spokes

Wheel Bag to accommodate four 24", 25" or 26" Wheels

