

STALLION

ADJUSTABLE FIXED FRAME DAY WHEELCHAIR



The stallion is a fixed frame wheelchair that has adjustable camber, multi-axle position and folding back. It is extremely versatile.

A sturdy lightweight wheelchair with all the adjustability you may require in a wheelchair.

This is a favourite everyday chair that can always be relied upon.

stallion
MELROSE LIFESTYLE WHEELCHAIR

Melrose
KIWI CONCEPT CHAIRS

STALLION

ADJUSTABLE FIXED FRAME DAY



Weight Limit: 115kg max, Custom weight option also available
Approx Weight: 6.3kg - Titanium standard 16" frame without wheels
7.0kg - Aluminium standard 16" frame without wheels

STANDARD SPECIFICATIONS

Aluminium Frame 90100

Adjustable Back Height and Angle

Folding or Fixed Back

Adjustable Foot-plate Height

Foot-plates: Swing-Away
Single Drop In
Fixed, Fixed Flip-Up

Adjustable Seat Rake, Camber and Balance Point

Castors: Roller-blade
4", 5", 6" Solid
4", 5", 6" Wide Solid

Quick release 24" Spoked Wheels on Alloy Rims

Tyres: 24 x 1 High Pressure Tread or Smooth
24 x 1 Solid

Brakes: Push to Lock Brakes (high or low mounted)

Alloy Scissor Brakes (high or low mounted)

Black Fabric Seat Sling, Back and Leg Strap or Heel Loops

Five Year Frame Warranty

STALLION BARIATRIC

Aluminium Frame 90100B

Seat Width: 18", 20", 22", 24", 26"

Wheels: Bariatric 24"

Weight Limit: 250kg

OPTIONS AND ACCESSORIES

Titanium Frame 90100T

Castors: 4", 5", 6" Solid Alloy
4", 5", 6" Soft Roll
6", 8" Pneumatic

Frog Leg Forks

Quad Release Axles

Wheels: Spoked 16", 18", 20", 22", 24", 25", 26"
BMX Wheels
Xcore Wheels
Spinergy Wheels
Soft Wheels
Loop Wheels

Push-rims: Rubber Coated - Red or Black
Vinyl
Foam
CarboLife
Natural Fit

Push handles: In Back Canes
Removable
Fold Down

Brakes: Composite Scissor Brakes

Eagle Scissor Brakes

Ds-Locks Brakes

Brake Extension Levers

Hill Climber Brakes

Anti-Tips with Plastic Wheels or Roller-blades

Spoke Guards: Clear or Black

Armrest: Tubular Swing Away

Desk or Full Length Pad

T Adjustable Height, with/without Skirt Guards

Adjustable Back Upholstery

Skirt Guards: Fabric Skirt Guards

Fabric covered ABS Skirt Guards

ABS Drop In Skirt Guards

Alloy or Carbon Fibre Skirt Guards

Under Chair Bag

Fabric or ABS Foot-plate Cover

Cushion, Impact Guards, Castor Covers