

TORPEDO

AFL SPORTS WHEELCHAIR



Built to cope with the big hits, these sports chairs have been designed to cope with the rigors of the fast expanding game of Aussie Rules. Developed with a robust frame built with reinforced bottom rails, this exciting new sport needs chairs with additional strength that keep going until the final whistle.

We've also added the highest quality components available,

many of these parts are built in-house by our team of qualified engineers using the latest CNC technology.

The Melrose Torpedo comes standard in lightweight aluminium, but can also be built in titanium for a super light super responsive frame.



TORPEDO

AFL SPORTS WHEELCHAIR

STANDARD SPECIFICATIONS

ALUMINUM FRAME

TITANIUM FRAME

Complete custom fit - built to individual specifications for size, rake and camber

Fixed camber bar between 12° and 20°

24", 25" or 26" spoked wheels with high pressure tyres on alloy rims

Hard anodised pushrims and quick release axles

Rollerblade castors

Adjustable height footplate

Cushion

Adjustable height back

Adjustable tension seat sling, back upholstery, fabric skirt guards, foot strap

OPTIONS AND ACCESSORIES

Tapered seat width and back

Adjustable tension back upholstery

Hip strap, back bar cover, cushions 1" 2" 3" or 4"

Double or single swivel anti-tips

Aluminium, Titanium or ABS skirt guards

Rubber coated Hand Rims - Red or Black

Spoke guards - clear or black

Tyres, smooth, high pressure IRC and Kenda Iron Cap

Super light weight Spinergy wheels - 24", 25", 26" or 700c with colour option carbon cable spokes

ABS or alloy footplate fixed or fully adjustable

Custom binding mounts

Various spoke guard designs:

- Polycarb
- Polyethylene
- Polyethylene / red rubber
- Alloy spoke guards

Wheel bag to accommodate (two or four) 24", 25" or 26" wheels

Bead blasted - titanium, or powder coated frame finish

Ergo Seating

