

Hornet Script Form

Fixed Frame Wheelchair 90460



CLIENT DETAILS	Date:
Name:	
Address:	
Phone / Mobile:	
Email:	

Funder:
Therapist:
Contact Ph:
Email:
Delivery Details:

Disability:	Trial Serial No.:
Weight: (actual) <input type="checkbox"/> <60kg <input type="checkbox"/> 60-100kg <input type="checkbox"/> >100kg	Height: Shoe Size: T-Shirt Size:

FABRICATION INFORMATION
Include details of any Power Adds and or front wheel/Freewheel (If necessary include diagrams and additional pages)

Melrose Consultant / Dealer:

Chair Serial No.:
Date Ordered:
Date Required:

Script checked and approved for build:	Date:
Name:	Signature:

Hornet Script Form

Fixed Frame Wheelchair

90460









Client Name: _____

Chair Serial No.: _____

MEASUREMENTS (Refer to diagrams)

Frame: Titanium 90460

A	Seat Width
B	Seat Length
	<input type="checkbox"/> Growth in Seat Length AD10700 (extends at rear)
	<input type="checkbox"/> 3/4 <input type="checkbox"/> 1-1/2 <input type="checkbox"/> 2-1/4 <input type="checkbox"/> 3
C1	Back Height
	<input type="checkbox"/> Fixed <input type="checkbox"/> Adjustable - Height
	Back Rest Style
	<input type="checkbox"/> Fixed <input type="checkbox"/> Folding <input type="checkbox"/> Folding/Locking AD11250
D	Back Angle°
E	Seat to Footplate
	<input type="checkbox"/> Fixed <input type="checkbox"/> Adjustable - Position
G	Frame Length
G1	Balance Point (COG)
H	Seat to Floor - Front
H1	Seat to Floor - Rear
I	Castor Pivot Centres
I1	Castor Pivot Position
J1	Footplate Width
K	Between Wheels
L	Camber
	<input type="checkbox"/> 0° <input type="checkbox"/> 1° <input type="checkbox"/> 2° <input type="checkbox"/> 3° <input type="checkbox"/> 4° <input type="checkbox"/> 5°
	Camber Bar Position
	<input type="checkbox"/> Fixed <input type="checkbox"/> Adjustable - Position
	Front Hanger Style
	<input type="checkbox"/> V Front <input type="checkbox"/> Pro T Front
	<input type="checkbox"/> Rolled Bottom Rail Option
	<input type="checkbox"/> Ergo Seat Option AD10220
M	Ergo Seat Length
M1	Ergo Back Height
M2	Ergo Back Angle°
	<input type="checkbox"/> Suspension Unit Option AD10200
	Pivots
	<input type="checkbox"/> STD <input type="checkbox"/> BULLET



Hornet Script Form

Fixed Frame Wheelchair

90460

Client Name: _____

Chair Serial No.: _____

FOOTPLATE	
<input type="checkbox"/> Standard	Std
<input type="checkbox"/> Angle adjustable	AD11520
<input type="checkbox"/> ABS on footplate	AD11500
<input type="checkbox"/> Fabric padded cover	AD11510
<input type="checkbox"/> Carbon Fibre cover	AD11506
<input type="checkbox"/> ABS with sides	AD11505
<input type="checkbox"/> Flip back	AD11452
<input type="checkbox"/> Custom:	

CASTORS	
<input type="checkbox"/> 3" Rollerblade	Std
<input type="checkbox"/> 4" x 1.25" wide solid composite	Std
<input type="checkbox"/> 5" x 1" solid composite	Std
<input type="checkbox"/> 6" x 1" solid composite	Std
<input type="checkbox"/> 3" x 1.5" alloy soft roll	AD11800
<input type="checkbox"/> 4" x 1.25" wide alloy	AD11710
<input type="checkbox"/> 4" x 1.5" alloy soft roll	AD11806
<input type="checkbox"/> 5" x 1" alloy	AD11740
<input type="checkbox"/> 5" x 1.5" wide solid composite	AD11720
<input type="checkbox"/> 5" x 1.5" wide alloy	AD11721
<input type="checkbox"/> 5" x 1.5" alloy soft roll	AD11810
<input type="checkbox"/> 6" x 1" alloy	AD11741
<input type="checkbox"/> 6" x 1.5" wide solid composite	AD11730
<input type="checkbox"/> 6" x 1.5" wide alloy	AD11742
<input type="checkbox"/> 6" x 1.5" alloy soft roll	AD11820
<input type="checkbox"/> 6" pneumatic	AD11750
<input type="checkbox"/> Custom:	

FORKS	
<input type="checkbox"/> Standard forks - black (for std pivots only)	Std
<input type="checkbox"/> Ultralight bearing forks (for bullet pivots only)	AD11910
<input type="checkbox"/> Frog Legs suspension	AD11950
<input type="checkbox"/> Frog Legs Unitine	AD11920
<input type="checkbox"/> Custom:	

BRAKES	
<input type="checkbox"/> Nil	
<input type="checkbox"/> Push to lock <input type="checkbox"/> Reverse handles	Std
<input type="checkbox"/> Scissor	Std
<input type="checkbox"/> Melrose L scissor	AD14522
<input type="checkbox"/> Melrose L scissor w/spring	AD14524
<input type="checkbox"/> Pull to lock	AD14510
<input type="checkbox"/> Hill climbers	AD14530
<input type="checkbox"/> Push to lock composite	AD14571
<input type="checkbox"/> Scissor lock composite	AD14570
<input type="checkbox"/> Eagle scissor <input type="checkbox"/> Black <input type="checkbox"/> Red <input type="checkbox"/> Blue	AD14573
<input type="checkbox"/> Superlock push to lock with brake extensions	AD14560
<input type="checkbox"/> D's locks	AD14600

BRAKES	
<input type="checkbox"/> Brake extensions <input type="checkbox"/> Left only <input type="checkbox"/> Right only	AD14550
<input type="checkbox"/> Custom:	

PUSH HANDLES	
<input type="checkbox"/> Nil	
<input type="checkbox"/> Height adj dual clamp	AD14120
<input type="checkbox"/> Foldable composite	AD14160
<input type="checkbox"/> Fold down alloy, with rubber grip	AD14161
<input type="checkbox"/> Fold down alloy, snaplock	AD14162
<input type="checkbox"/> In back canes	AD14140
<input type="checkbox"/> Removable on grab bar	AD14100
<input type="checkbox"/> Ratchet Stroller	AD14300
PUSH HANDLE(S) SET AT HEIGHT:	
<input type="checkbox"/> Custom:	

ANTI TIPS	
<input type="checkbox"/> Nil	
<input type="checkbox"/> Plastic wheels	AD15000
<input type="checkbox"/> Stainless steel with plastic wheels	AD15025
<input type="checkbox"/> Rollerblades	AD15030
<input type="checkbox"/> Stainless steel with rollerblades	AD15026
<input type="checkbox"/> Latch flip up	AD15022
<input type="checkbox"/> Single centre mount with plastic wheel	AD15070
<input type="checkbox"/> Single centre mount with roller blade	AD15080
<input type="checkbox"/> Double mount one piece swivel anti tip with rollerblade	AD15150
<input type="checkbox"/> Swing under plastic wheels	AD15060
<input type="checkbox"/> Custom:	

ARMRESTS	
<input type="checkbox"/> Nil	
<input type="checkbox"/> T height adjustable desk pad armrest	AD15260
<input type="checkbox"/> T height adjustable full length pad armrest	AD15261
<input type="checkbox"/> T height adjustable desk pad armrest with ABS skirt guard	AD15280
<input type="checkbox"/> T height adjustable full length pad armrest with ABS skirt guard	AD15281
<input type="checkbox"/> Swing away, foam tube	AD15200
<input type="checkbox"/> Swing away, desk pad	AD15210
<input type="checkbox"/> Swing away, full pad	AD15220
<input type="checkbox"/> Flip back, foam tube on back cane	AD15225
ARMRESTS SET AT HEIGHT:	
<input type="checkbox"/> Custom:	

Hornet Script Form

Fixed Frame Wheelchair 90460

Client Name: _____

Chair Serial No.: _____

SKIRT GUARDS

NOTE: Folding Skirt Guards Limit Brake Options *

<input type="checkbox"/>	Nil	
<input type="checkbox"/>	Fabric skirt guards	AD16100
<input type="checkbox"/>	ABS fabric covered	AD16120
<input type="checkbox"/>	ABS drop in	AD16110
<input type="checkbox"/>	Carbon fibre drop in	AD16162
<input type="checkbox"/>	Alloy fixed	AD16130
<input type="checkbox"/>	Alloy folding	AD16140
<input type="checkbox"/>	Carbon fibre fixed	AD16161
<input type="checkbox"/>	Carbon fibre folding with back *	AD16160
<input type="checkbox"/>	Carbon fibre slotted folding (through)	AD16158
<input type="checkbox"/>	Carbon fibre fold in (on to sling)	AD16172
<input type="checkbox"/>	Carbon fibre flip back	AD16174
<input type="checkbox"/>	Alloy flip back	AD16141
<input type="checkbox"/>	Carbon fibre drop in fenders	AD16163
<input type="checkbox"/>	Alloy drop in mud guard	AD16164
SKIRT GUARDS SET AT HEIGHT:		
<input type="checkbox"/>	Custom:	

UPHOLSTERY - Seating

<input type="checkbox"/>	Sling - screwed down	Std
<input type="checkbox"/>	Sling - full straps	Std
Sling Velcro: <input type="checkbox"/> Hook <input type="checkbox"/> Loop		
<input type="checkbox"/>	Solid alloy base	AD18030
<input type="checkbox"/>	Solid carbon base	AD18031
<input type="checkbox"/>	Tension adjustable back	AD16010
<input type="checkbox"/>	Standard back	Std
<input type="checkbox"/>	No back upholstery	
<input type="checkbox"/>	To fit solid back:	
<input type="checkbox"/>	Cushion 1"	AD16201
<input type="checkbox"/>	Cushion 2"	AD16202
<input type="checkbox"/>	Cushion 3"	AD16203
<input type="checkbox"/>	Custom:	

UPHOLSTERY - Bags

<input type="checkbox"/>	Bag on front hanger - left	AD16620
<input type="checkbox"/>	Bag on front hanger - right	AD16620
<input type="checkbox"/>	Bag on back of chair	AD16680
<input type="checkbox"/>	Bag under chair	AD16600
<input type="checkbox"/>	Bag under chair - Net	AD16650
<input type="checkbox"/>	Bag under chair - quick release, small	AD16601
<input type="checkbox"/>	Wallet / phone holder	AD16622
<input type="checkbox"/>	Drink bottle holder	AD16624
<input type="checkbox"/>	Custom:	

UPHOLSTERY - Leg Supports

<input type="checkbox"/>	Leg strap	Std
<input type="checkbox"/>	Padded calf strap	AD16320
<input type="checkbox"/>	Body Point calf strap - Size:	

UPHOLSTERY - Impact Guards

(Impact guards made only if suitable after fitting anti tips, push handles etc.)

<input type="checkbox"/>	Front hanger impact guards - Webbing Fabric	AD16400
<input type="checkbox"/>	Front hanger impact guards - Neoprene	AD16401
<input type="checkbox"/>	Front hanger impact guards - Leather	AD16402
<input type="checkbox"/>	Grab bar impact guards	AD16420
<input type="checkbox"/>	Bottom bar impact guards	AD16422

SEAT BELT

<input type="checkbox"/>	Car buckle	AD16522
<input type="checkbox"/>	Plastic buckle	AD16500
<input type="checkbox"/>	Body Point - Size/Style:	

ACCESSORIES

<input type="checkbox"/>	Tray - pronged	AD17010
<input type="checkbox"/>	Tray - slides over armrests	AD17000
<input type="checkbox"/>	Tie down mounts	AD10320

FRAME COLOURS

Beadblast
StdPolished
AD10120

(Powder coat colours are not recommended on Titanium frames.)

Hornet Script Form

Fixed Frame Wheelchair 90460

Client Name: _____

Chair Serial No.: _____

STANDARD REAR WHEELS (for BMX - see Spare Wheels)

		22	24	25	26
Spoked	Std	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

PERFORMANCE REAR WHEELS

		22	24	25	26
Spinergy L Spox (18 spoke)	AD13300		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Spinergy LX Lite (12 spoke)	AD13310	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Spinergy LX with FlexRim (black only)	AD13315		<input type="checkbox"/>	<input type="checkbox"/>	
Spinergy XLX (12 spoke cross laced)	AD13312		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Spinergy Blade (18 spoke)	AD13313	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Spinergy Carbon CLX (18 spoke)	AD13318		<input type="checkbox"/>	<input type="checkbox"/>	
Spinergy Sport SLX (24 spoke large hub heavy duty - black only)	AD13420		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Spinergy Sport XSLX (24 spoke cross laced - black only)	AD13314		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Custom wheel (specify below)		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
-					

SPINERGY SPOKE COLOUR (Optional)

(Colour spokes not available on Carbon, Sports and FlexRim wheels)

Spinergy L Spox	<input type="checkbox"/> Black <input type="checkbox"/> Red <input type="checkbox"/> Blue <input type="checkbox"/> White <input type="checkbox"/> Yellow
Spinergy LX, Spinergy Blade	<input type="checkbox"/> Black <input type="checkbox"/> Red <input type="checkbox"/> Blue <input type="checkbox"/> White <input type="checkbox"/> Yellow
	<input type="checkbox"/> Green <input type="checkbox"/> Purple <input type="checkbox"/> Pink <input type="checkbox"/> Orange
Custom -	

AXLES

<input type="checkbox"/> Quick Release	Std	
<input type="checkbox"/> Quad Release	AD12500	

TYRES

		22	24	25	26
1" HP Smooth	Std	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
1" HP Tread	Std		<input type="checkbox"/>		
1" Solid	Std	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
1-3/8" Tread	Std	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>
1-3/8" Solid	Std	<input type="checkbox"/>	<input type="checkbox"/>		
1-3/8" Kenda HP Tread	AD13710		<input type="checkbox"/>		
1-3/8" Knobbly V-Trak	AD13715		<input type="checkbox"/>		
1" Kenda Ironcap Blue	AD13700		<input type="checkbox"/>	<input type="checkbox"/>	
1" Kenda Konstrictor	AD13720		<input type="checkbox"/>	<input type="checkbox"/>	
1" Schwalbe Right Run	AD13740	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
1" Schwalbe Marathon Plus Evo	AD13750		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
1-3/8" Schwalbe Marathon Plus Evo	AD13751		<input type="checkbox"/>		
Custom -		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

PUSHRIMS

		22	24	25	26
<input type="checkbox"/> Nil					
Hard Anodised	Std	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Hard Anodised 1"	AD13855		<input type="checkbox"/>		
Vinyl	AD13830	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Vinyl 1"	AD13859		<input type="checkbox"/>	<input type="checkbox"/>	
Foam	AD13840	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Foam 1"	AD13841		<input type="checkbox"/>		
Rubber coated (black)	AD13850	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Rubber coated 1"	AD13857		<input type="checkbox"/>	<input type="checkbox"/>	
Surge Lite	AD13864	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Natural Fit Lite	AD13860	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/> With thumb grip <input type="checkbox"/> Without thumb grip					
Max Grippe	AD13870		<input type="checkbox"/>	<input type="checkbox"/>	
Stainless Steel	AD13810	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Spinergy Stainless Textured	AD13811	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Spinergy Max Stainless Textured (Ergo)	AD13812		<input type="checkbox"/>	<input type="checkbox"/>	
Q Grips	AD13865		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Titanium	AD13820		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Carbolife	AD13868		<input type="checkbox"/>	<input type="checkbox"/>	
<input type="checkbox"/> Curve <input type="checkbox"/> Curve Large					
<input type="checkbox"/> Curve Tetra <input type="checkbox"/> Curve Tetra Large					
<input type="checkbox"/> Gekko <input type="checkbox"/> Gekko S					
<input type="checkbox"/> Quadro					

PUSHRIMS - Position

<input type="checkbox"/> Std	<input type="checkbox"/> In Close	<input type="checkbox"/> On to rim
------------------------------	-----------------------------------	------------------------------------

SPOKEGUARDS

<input type="checkbox"/> Clear polycarb	AD13900	
<input type="checkbox"/> Black ABS	AD13901	
<input type="checkbox"/> Custom design (specify below)	AD13910	
Design name -		

SPARE REAR WHEELS

<input type="checkbox"/> BMX Wheels	<input type="checkbox"/> Extra Wheels
Description:	
Size:	
Tyre:	
Pushrim:	
Setting:	
To interchange with:	