

CLIENT DETAILS		Date:		Funder:		
Name:				Therapist:		
Address:				Contact Ph:		
				Email:		
				Delivery Details:		
Phone / Mobile:						
Email:						
Disability:				Trial Serial No.:		
Weight: (actual)	☐ <60kg	☐ 60-100kg	☐ >100kg	Height:	Shoe Size:	T-Shirt Size:
FABRICATION INFORM				Melrose Consultant / D	Dealer:	
nclude details of any Power Add: If necessary include diagrams and addition	s and or front w	/heel/Freewheel				

Script checked and approved for build:

Date:

Signature:

Chair Serial No.:


Date Ordered:

Date Required:



MEASUREMENTS (Refer to diagrams)						
Frame: Titanium 90356TS Aluminium 90356S						
Α	Seat Width 15.5" — 20"					
В	Seat Length					
	Growth in Seat Length AD10700 (extends at rear)					
	☐ 3/4" ☐ 1-1/2" ☐ 2-1/4"					
C1	Back Height					
	☐ Fixed ☐ Adjustable - Height					
	Back Rest Style					
	☐ Fixed ☐ Folding					
D	Back Angle°					
Е	Seat to Footplate					
	☐ Fixed ☐ Adjustable - Position					
G	Frame Length					
G1	Balance Point (COG)					
Н	Seat to Floor - Front					
H1	Seat to Floor - Rear Range 16" - 18.5"					
	To be confirmed on CAD drawing with tyre choice and setup.					
I	Castor Pivot Centres					
I1	Castor Pivot Position					
	Standard Position Pivot Forward Position					
J1	Footplate Width					
L	Camber					
	□ 0° □ 1° □ 2° □ 3° □ 4° □ 5°					
	Front Hanger Style					
	☐ V Front ☐ Pro T Front					
М	Ergo Seat Length					
	Pivots					

☐ Standard ☐ Bullet





FOOTPLATE	
Standard	Std
Angle adjustable	AD11520
ABS on footplate	AD11500
Fabric padded cover	AD11510
Carbon Fibre cover	AD11506
ABS with sides	AD11505
Custom:	
CASTORS	
4" x 1.5" alloy soft roll	AD11806
5" x 1.5" wide solid composite	AD11720
5" x 1.5" wide alloy	AD11721
5" x 1.5" alloy soft roll	AD11810
6" x 1.5" wide solid composite	AD11730
6" x 1.5" wide alloy	AD11742
6" x 1.5" alloy soft roll	AD11820
6" pneumatic	AD11750
☐ 175 x 50 offroad	AD11784
	'
FORKS	
Standard forks - black (for std pivots only)	Std
Ultralight bearing forks (for bullet pivots only)	AD11910
Frog Legs suspension	AD11950
Hercules forks (for bullet pivots only)	AD11952
Custom:	
BRAKES	
Nil	
Push to lock Reverse handles	Std
Scissor	Std
Pull to lock	AD14510
Hill climbers	AD14530
Push to lock composite	AD14571
Scissor lock composite	AD14570
☐ Eagle scissor ☐ Black ☐ Red ☐ Blue	AD14573
☐ Brake extensions ☐ Left only ☐ Right only	AD14550
Custom:	

Chair	Serial No.:	
PU	SH HANDLES	
	Nil	
	Height adj dual clamp	AD14120
	Foldable composite	AD14160
	Fold down alloy, with rubber grip	AD14161
	Fold down alloy, snaplock	AD14162
	In back canes	AD14140
	Removable on grab bar	AD14100
	PUSH HANDLE(S) SET AT HEIGHT:	
	Custom:	

Client Name:

Client Name:



		Chair Serial No.:
SKIRTGUARDS		UPHOLSTERY - Leg Su
NOTE: Folding Skirt Guards Limit Brake Options *		Leg strap
□ Nil		Padded calf strap
Fabric skirt guards	AD16100	Body Point calf strap - Size
☐ ABS fabric covered	AD16120	
☐ ABS drop in	AD16110	UPHOLSTERY - Impac
Carbon fibre drop in	AD16162	(Impact guards made only if su
☐ Alloy fixed	AD16130	
☐ Alloy folding	AD16140	Front hanger impact guard
☐ Carbon fibre fixed	AD16161	Front hanger impact guard
☐ Carbon fibre slotted folding (through)	AD16158	Front hanger impact guard
☐ Carbon fibre drop in fenders	AD16163	Grab bar impact guards
Alloy drop in mud guard	AD16164	Bottom bar impact guards
	·	
SKIRT GUARDS SET AT HEIGHT:		SEAT BELT
Custom:		Cambrodia
		Car buckle
UPHOLSTERY - Seating		Plastic buckle
		Body Point - Size/Style:
Sling - screwed down	Std	
Sling - full straps	Std	
Sling Velcro: Hook Loop		
Solid alloy base	AD18030	
Solid carbon base	AD18031	
Tension adjustable back	AD16010	
Standard back	Std	
No back upholstery		
To fit solid back:		
Cushion 1"	AD16201	
Cushion 2"	AD16202	
Cushion 3"	AD16203	
Custom:		
UPHOLSTERY - Bags		
Bag on front hanger - left	AD16620	
Bag on front hanger - right	AD16620	
☐ Bag on back of chair	AD16680	
Bag under chair	AD16600	
Bag under chair - Net	AD16650	
Bag under chair - quick release, small	AD16601	
Wallet / phone holder	AD16622	
☐ Drink bottle holder	AD16624	

Custom:

Chair	Serial No.:				
UPHOLSTERY - Leg Supports					
	Leg strap	Std			
	Padded calf strap	AD16320			
	Body Point calf strap - Size:				
UP	HOLSTERY - Impact Guards				
(Impact guards made only if suitable after fitting anti tips, push handles etc.)					
	Front hanger impact guards - Webbing Fabric	AD16400			
	Front hanger impact guards - Neoprene	AD16401			
	Front hanger impact guards - Leather	AD16402			
	Grab bar impact guards	AD16420			
	Bottom bar impact guards	AD16422			
SEAT BELT					
	Car buckle	AD16522			
	Plastic buckle	AD16500			
	Body Point - Size/Style:				



22 24 25

					Chair Serial No.:		
STANDARD REAR WHEELS (f	or BMX - so	ee Spare	Whe	els)	PUSHRIMS		
(1		•		-			
		22	24	25			507
Spoked	Std				□ Nil		
					Hard Anodised	Std	
PERFORMANCE REAR WHEI	ELS				Hard Anodised 1"	AD13855	
					Rubber coated (black)	AD13850	
		22	24	25	Rubber coated 1"	AD13857	
Spinergy LX Lite (12 spoke)	AD13310		ᆜ	ᆜ	Surge Lite	AD13864	
Spinergy XLX (12 spoke cross laced)	AD13312		ᆜ	ᆜ	Natural Fit Lite	AD13860	
Spinergy Blade (18 spoke)	AD13313		Ц	Щ	☐ With thumb grip ☐ Without thun		
Spinergy Carbon CLX (18 spoke)	AD13318			ᆜ	Stainless Steel	AD13810	
Spinergy Sport SLX (24 spoke large	AD13420		Ш	Ш	Spinergy Stainless Textured	AD13811	
hub heavy duty - black only)	1012211				Spinergy Max Stainless Textured (Ergo)	AD13812	
Spinergy Sport XSLX (24 spoke cross laced - black only)	AD13314			Ш	Q Grips	AD13865	
Spinergy BMX (507 - 18 spoke	AD13339				Titanium	AD13820	
- black only)	ADISSSS		Ш		Carbolife	AD13868	
Custom wheel (specify below)			П		Curve Curve Large		
Custom wheel (specify below)			Ш	ш	☐ Curve Tetra ☐ Curve Tetra Large ☐ Gekko ☐ Gekko S		
					Quadro		
		-					
SPINERGY SPOKE COLOUR	(Optional)					
(Colour spokes not available on Carbor	n, Sports and F	lexRim whe	els)		PUSHRIMS - Position		
☐ Std ☐ In Close ☐ On to rim							
Spinergy L Spox							
	te Yellow				SPOKEGUARDS		
Spinergy LX, Spinergy Blade	ta □Valla				SPOREGUARDS		
	te Yellow				☐ Clear polycarb	AD13900	
	ige				☐ Black ABS	AD13901	
Custom -					Custom design (specify below)	AD13910	
Design name -							
AXLES							
Quick Release	Std				SPARE REAR WHEELS		
Quad Release	AD12500				☐ BMX Wheels ☐ E.	xtra Wheels	
					Description:		
TYRES					bescription.		
TINES					Size:		
		22	24	25	3.20.		
1.95" BMX Tread (507)	Std				Tyre:		
2.125" BMX Tread (507)	Std				,,		
1-3/8" Knobbly V-Trak	Std				Pushrim:		
1" Kenda Ironcap Blue	Std						
2.1" Kenda Nevegal	Std				Setting:		
1" Kenda Konstrictor	AD13720				, 3		
					To interchange with:		
1" Schwalbe Marathon Plus Evo	AD13750				3 '		
1-3/8" Schwalbe Marathon Plus Evo	AD13751						
2" Schwalbe Billy Bonker (507)	AD13752						

Client Name:

2.1" Schwalbe Smart Sam (507)


AD13753




Client Name:		
Chair Serial No.:		

FRAME COLOURS





GLASS COLOURS




CANDY COLOURS



CANDY COLOURS (SPECIAL)





THIS IS A GUIDE TO COLOURS ONLY, HOW COLOURS APPEAR ON YOUR MONITOR OR WHEN PRINTED OUT MAY VARY.

Talk to us if you want your Titanium frame coloured or polished. Colour availability is subject to change without notice.