

# Stallion Script Form

Fixed Frame Adjustable Wheelchair 90100



**To Order** Fill in the frame dimensions and features required. For variations describe in detail on a separate page. If you have any further questions contact us at the address at the bottom of the page.

**SEND FORM**  
If online, click here to send via e-mail

Date \_\_\_\_\_

Name \_\_\_\_\_

Address \_\_\_\_\_

Phone \_\_\_\_\_

Fax / Cell \_\_\_\_\_

Email \_\_\_\_\_

Delivery details \_\_\_\_\_

Therapist details \_\_\_\_\_

Required by \_\_\_\_\_

Disability \_\_\_\_\_ Weight \_\_\_\_\_ Height \_\_\_\_\_ Shoe Size \_\_\_\_\_ T-Shirt Size \_\_\_\_\_

## MEASUREMENTS (Refer to diagram)

**Materials**  Aluminium (std)  Titanium **90100T**

**A** Seat width 14"  15"  16"  17"   
18"  19"  20"

**A** Custom seat width **AD10100**

**B** Seat depth 14"  16"  18"  20"

**B** Custom seat depth **AD10100**

**B** Growth in seat depth **AD10900**

**C** Back height \_\_\_\_\_

**D** Back angle° \_\_\_\_\_

**E** Seat to footplate \_\_\_\_\_

**F** Seat to floor - front \_\_\_\_\_

**G** Seat to floor - rear \_\_\_\_\_

**H** Balance Point \_\_\_\_\_

**I** Castor Pivot Position \_\_\_\_\_

**J** Castor Pivot Centres \_\_\_\_\_

**Max User Weight**  
115kg for up to 17" seat width  
165kg for 18" seat width and above


**K** Camber (0°-5°) \_\_\_\_\_

**L** Hanger Angle  
 70°  80° (std)  90°

**Back Style**  
 Folding Back  Fixed Back  
 **Dynamic Rocker Back** **AD10330**  
Additional user weight up to 165kg  
 Reinforced aluminium frame **AD10400**  
 Reinforced titanium frame **AD10410**

## ADDITIONAL NOTES

(you can type here)



# Stallion Script Form

Fixed Frame Adjustable Wheelchair


Name \_\_\_\_\_

**HANGER TYPE**

Swing Away



Fixed




Drop In



**FOOT-PLATE WIDTH**

Width for Pro T and V style foot-plates




When measuring across feet you must add 2" (50mm)

Front taper standard foot-plate width measurements


SEAT WIDTH	14"	16"	18"	20"
FOOT-PLATE WIDTH	12"	12"	13"	15"
	(305mm)	(305mm)	(340mm)	(381mm)

**SWING AWAY FOOTREST**


Straight Swing Away




Straight Swing Away



V Swing Away



V Swing Away



Straight Hangers

2 piece, flip up angle adjustable alloy foot-plates (std)

V Hangers

1 piece, flip up angle adjustable alloy foot-plate, opens to the left **AD11555**

1 piece, flip up angle adjustable alloy foot-plate, opens to the right (pictured) **AD11555**

**FIXED FOOTREST**

Straight Fixed  
D style  
foot-plate



Straight Fixed  
2 flip up



Straight Fixed  
1 flip up



Pro T Fixed  
D style  
foot-plate



Pro T Fixed  
2 flip up



Pro T Fixed  
1 flip up



V Fixed  
D style  
foot-plate



V Fixed  
2 flip up



V Fixed  
1 flip up



Straight Hangers

D style foot-plate (std)

ABS on foot-plate (std)

Pro T Hangers

Angle adjustable foot-plate (std)

V Hangers

2 piece, flip up angle adjustable alloy foot-plates (std)

1 piece, flip up angle adjustable alloy foot-plate, opens to the left **AD11550**

1 piece, flip up angle adjustable alloy foot-plate, opens to the right (pictured) **AD11550**

**DROP IN FOOTREST**

Straight Drop In  
D style  
foot-plate



Pro T Drop In  
D style  
foot-plate



Straight Hangers

D style foot-plate (std)

ABS on foot-plate (std)

Pro T Hangers

Angle adjustable foot-plate (std)

**ELEVATING LEG REST**

Made for the Swing Away front only

Left **AD11610**

Right **AD11610**

Both **AD11600**

**STUMP SUPPORT**

Made for the Fixed and Swing Away fronts only

Left **AD11650**

Right **AD11650**

Both **AD11651**

Fixed

Removable

# Stallion Script Form

## Fixed Frame Adjustable Wheelchair

Name \_\_\_\_\_

**CASTORS**

- |  |  |  |   |
|--|--|--|---|
| <input type="checkbox"/> 4" solid composite (std)        | <input type="checkbox"/> 5" solid composite (std)          | <input type="checkbox"/> 6" solid composite (std)          | <input type="checkbox"/> 8" solid composite AD11760     |
| <input type="checkbox"/> 4" wide solid composite (std)   | <input type="checkbox"/> 5" solid alloy AD11740            | <input type="checkbox"/> 6" solid alloy AD11742            | <input type="checkbox"/> 8" floater composite AD11761   |
| <input type="checkbox"/> 4" x 1.25" solid alloy AD11710  | <input type="checkbox"/> 5" x 1.5" solid composite AD11720 | <input type="checkbox"/> 6" x 1.5" solid composite AD11730 | <input type="checkbox"/> 8" pneumatic composite AD11780 |
| <input type="checkbox"/> 4" wide soft roll alloy AD11806 | <input type="checkbox"/> 5" x 1.5" solid alloy AD11740     | <input type="checkbox"/> 6" x 1.5" solid alloy AD11742     |   |
|  | <input type="checkbox"/> 5" wide soft roll alloy AD11810   | <input type="checkbox"/> 6" wide soft roll alloy AD11820   |   |
|  |  | <input type="checkbox"/> 6" pneumatic composite AD11750    |   |

**FORKS**

- 
- Standard forks (std)
- 
- Other \_\_\_\_\_

**PUSH HANDLES**

- |                              |  |   |   |
|------------------------------|--|---|---|
| <input type="checkbox"/> Nil | <input type="checkbox"/> Removable on grab bar AD14100 | <input type="checkbox"/> In back canes AD14140          | <input type="checkbox"/> Stroller D AD14200                 |
|                              | <input type="checkbox"/> Dual clamp adj height AD14120 | <input type="checkbox"/> Foldable in back canes AD14160 | <input type="checkbox"/> Stroller handles removable AD14240 |

**BRAKES**

- |   |   |   |  |
|---|---|---|--|
| <input type="checkbox"/> Nil                                | <input type="checkbox"/> Pull to lock AD14510     | <input type="checkbox"/> Reverse handles                | <input type="checkbox"/> Footbrake AD14660 |
| <input type="checkbox"/> Push to lock (std)                 | <input type="checkbox"/> Brake extensions AD14550 | <input type="checkbox"/> Scissor AD14520                |  |
| <input type="checkbox"/> Push to lock (long handle) AD14551 | <input type="checkbox"/> Hill climbers AD14530    | <input type="checkbox"/> Scissor lock composite AD14570 |  |

**ANTI-TIPS**

- |  |  |  |  |
|--|--|--|--|
| <input type="checkbox"/> Nil               | <input type="checkbox"/> With plastic wheels (std) | <input type="checkbox"/> Stainless Steel with plastic wheels AD15025   | <input type="checkbox"/> Adjustable anti tips AD15010  |
| <input type="checkbox"/> Kick tube AD15040 | <input type="checkbox"/> With rollerblades AD15030 | <input type="checkbox"/> Flared w/plastic wheels AD15020<br><small>(only for w/chairs smaller than "16 wide)</small> | <input type="checkbox"/> Swing under anti tips AD15060 |

**ARMRESTS**

- |   |  |
|---|--|
| <input type="checkbox"/> Nil                                    | <input type="checkbox"/> T adjustable height with desk pad armrest AD15260                 |
| <input type="checkbox"/> Foam tube swing away, lift out AD15200 | <input type="checkbox"/> T adjustable height with full length pad armrest AD15261          |
| <input type="checkbox"/> Desk pad swing away, lift out AD15210  | <input type="checkbox"/> T adjustable height desk pad armrest with ABS skirt guard AD15280 |
| <input type="checkbox"/> Full pad swing away, lift out AD15220  | <input type="checkbox"/> T adjustable height full pad armrest with ABS skirt guard AD15281 |

**UPHOLSTERY**

- |  |  |  |  |
|--|--|--|--|
| <input type="checkbox"/> Sling (std)                   | <input type="checkbox"/> To fit solid back _____ | <input type="checkbox"/> Cushion 2" AD16202      | <input type="checkbox"/> Bag on back of chair AD16680      |
| <input type="checkbox"/> No back upholstery            | _____  | <input type="checkbox"/> Cushion 3" AD16203      | <input type="checkbox"/> Bag on front left hanger AD16620  |
| <input type="checkbox"/> Tension adjustable back (std) | <input type="checkbox"/> Cushion 1" AD16201      | <input type="checkbox"/> Bag under chair AD16600 | <input type="checkbox"/> Bag on front right hanger AD16620 |

**SKIRT GUARDS**

- |  |   |  |  |
|--|---|--|--|
| <input type="checkbox"/> Nil                         | <input type="checkbox"/> ABS drop in AD16110        | <input type="checkbox"/> Alloy with fixed back AD16130   | <input type="checkbox"/> Carbon Fibre folding back AD16160 |
| <input type="checkbox"/> Fabric skirt guards AD16100 | <input type="checkbox"/> ABS fabric covered AD16120 | <input type="checkbox"/> Alloy with folding back AD16140 |  |

**STRAPS / BELTS, GUARDS, GLOVES**

- |  |  |   |
|--|--|---|
| <input type="checkbox"/> Leg strap (std)                             | <input type="checkbox"/> Impact guards hangers (std)     | <input type="checkbox"/> Seat belt with plastic buckle AD16500                |
| <input type="checkbox"/> Heel loops (std)                            | <input type="checkbox"/> Impact guards back bars AD16420 | <input type="checkbox"/> Seat belt (adult) with metal car type buckle AD16522 |
|  | <input type="checkbox"/> Castor pivot covers AD16410     | <input type="checkbox"/> Hip strap 4 point belt AD16540                       |
| <input type="checkbox"/> Gloves<br><small>(size &amp; style)</small> | <input type="checkbox"/> Fabric footplate cover AD11510  |   |

**TRAY**

- 
- With prong receivers AD17010
- 
- 
- Slides over armrests AD17000

**TIE DOWN MOUNTS**

- 
- 2 front and 2 rear AD10320

# Stallion Script Form

## Fixed Frame Adjustable Wheelchair

Name \_\_\_\_\_

**WHEELS****16" - 22"**

- 16" x 1.75" Composite Wheels  
 18" x 1" Spoked Wheels  
 18" x 1.75" Composite Wheels  
 20" x 1" Spoked Wheels  
 20" x 1-3/8" Spoked Wheels  
 20" x 1.75" Spoked Wheels  
 22" x 1" Spoked Wheels  
 22" x 1-3/8" Spoked Wheels  
 22" Spinerger LX Wheels (12 spoke - black only) **AD13310**

**24"**

- Spoked Wheels (std)  
 Spinerger Day Wheels (18 spoke) **AD13300** (Choose your spoke colour) **Bk R Bl W Y**  
 Spinerger LX Wheels (12 spoke) **AD13310** (Choose your spoke colour) **Bk R Bl W Y Gr Pk O**  
 Spinerger LX Wheels with FlexRim (12 spoke - black only) **AD13315**  
 X Core Wheels (5 spoke) **AD13108**  
 Dino Wheels **AD13360** (Choose your spoke colour) **Bk R Bl Sv Go Gr Pk O LB Pu**  
 24" x 1.75" BMX Wheels (std)

**25"**

- Spoked Wheels (std)  
 Spinerger Day (18 spoke) **AD13300** (Choose your spoke colour) **Bk R Bl W Y**  
 Spinerger LX (12 spoke) **AD13310** (Choose your spoke colour) **Bk R Bl W Y Gr Pk O**  
 Spinerger LX Wheels with FlexRim (12 spoke - black only) **AD13315**  
 X Core Wheels (3 spoke) **AD13200**  
 Dino Wheels **AD13360** (Choose your spoke colour) **Bk R Bl Sv Go Gr Pk O LB Pu**  
 24" x 1.95" BMX Wheels (std)

**26"**

- Spoked Wheels (std)  
 Spinerger Day (18 spoke) **AD13300** (Choose your spoke colour) **Bk R Bl W Y**  
 Spinerger LX (12 spoke) **AD13310** (Choose your spoke colour) **Bk R Bl W Y Gr Pk O**  
 26" x 1.75" BMX Wheels (std)

 Spare Wheel(s) Type: \_\_\_\_\_

Tyres: \_\_\_\_\_

Pushrims: \_\_\_\_\_

Size: \_\_\_\_\_

**TYRES****16" - 22"**

- 16" x 1.75" Tread (std)  
 18" x 1" HP Smooth (std)  
 18" x 1.75" Tread (std)  
 20" x 1.75" Tread (std)  
 20" x 1" HP Smooth (std)  
 20" x 1-3/8" Tread (std)  
 20" x 1" Solid (std)  
 22" x 1" HP Smooth (std)  
 22" x 1-3/8" Tread (std)  
 22" x 1" Solid (std)

**24"**

- 24" x 1" HP Smooth (std)  
 24" x 1" HP Tread (std)  
 24" x 1" Solid (std)  
 24" x 1-3/8" Tread (std)  
 24" x 1-3/8" Solid (std)  
 24" x 1-3/8" Kenda HP Tread **AD13710**  
 24" x 1-3/8" Primo V-Trak Knobbly **AD13715**  
 24" x 1" Kenda Ironcap Blue **AD13700**  
 24" x 1" Kenda Konstrictor **AD13720**  
 24" x 1" Schwalbe Right Run **AD13740**  
 24" x 1" Schwalbe Marathon Evo **AD13750**  
 24" x 1.75" BMX Tread (std)

**25"**

- 25" x 1" HP Smooth (std)  
 25" x 1" Solid (std)  
 25" x 1" Kenda Ironcap Blue **AD13700**  
 25" x 1" Kenda Konstrictor **AD13720**  
 25" x 1" Schwalbe Right Run **AD13740**  
 25" x 1" Schwalbe Marathon Evo **AD13750**  
 24" x 1.95" BMX Tread (std)

**26"**

- 26" x 1" HP Smooth (std)  
 26" x 1-3/8" Tread (std)  
 26" x 1" Solid (std)  
 26" x 1" Schwalbe Right Run **AD13740**  
 26" x 1" Schwalbe Marathon Evo **AD13750**  
 26" x 1.75" BMX Tread (std)

**AXLES**

- Quick release axles (std)       Quad release **AD12500**

**SPOKE GUARDS**

- Clear polycarb **AD13900**       Black ABS **AD13901**

**PUSHRIMS**

- Anodised (std)       Surge **AD13864**       Natural Fit no thumb grips **AD13860**       Spinerger FlexRim see Wheels  
 Rubber Coated (black) **AD13850**       Max Grippe **AD13870**       Natural Fit with thumb grips **AD13860**  
 Vinyl **AD13830**       Capstan No. of projections: \_\_\_\_\_ Angle: \_\_\_\_\_ Length: \_\_\_\_\_

**PUSHRIMS - POSITION**

- Standard       In close       Onto rim



Name \_\_\_\_\_

### ATTENTION

**THIS IS A GUIDE TO COLOURS ONLY, how colours appear on your monitor or when printed out may vary.**

Talk to us if you want your Titanium frame coloured or polished. *Colour availability is subject to change without notice.*

#### STANDARD & ENABLE CLIENT COLOURS



#### GLASS COLOURS



#### CANDY COLOURS



#### CANDY COLOURS (SPECIAL)



#### FLUORESCENT COLOURS

